

# Multicolor Scan Laser Photocoagulator MC-500 Vixi Multicolor Laser Photocoagulator MC-500



THE ART OF EYE CARE



# The Versatile Laser Photocoagulator




## Multicolor on Modular Architecture

#### Multicolor Laser for Multiple Applications

The MC-500 Vixi / MC-500 enables efficient photocoagulation even through opaque media. In cases of cataract, better penetration is achieved with the yellow (577 nm) laser compared to the green (532 nm) laser.

In eyes with retinal hemorrhage, better penetration is achieved with the red (647 nm) laser.



#### ● 532 nm


The 532 nm (green) is the most common wavelength for treating retinal pathology.

#### 577 nm

The 577 nm (yellow) laser is minimally absorbed by xanthophyll and is well absorbed by oxygenated hemoglobin compared to 561 nm and 568 nm lasers making it the wavelength of choice for lesions close to the macula. This wavelength has plentiful results achieved with the dye lasers.

#### 647 nm

The 647 nm (red) wavelength has been historically used in krypton lasers. This wavelength is used for photocoagulation of deep choroidal pathology.



Reference: Folia Ophthalmol. Jpn. 40(5)1128-1133, 1989




#### Selectable Laser Color Configuration

The MC-500 Vixi / MC-500, with its user friendly design, allows the selection of one, two, or three wavelengths, among green, yellow, and red. It enables the freedom to select the necessary color or combination of colors to increase efficiency of treatment.

| Three-color selection |  |   |          |  |
|-----------------------|--|---|----------|--|
| Two-color selection   |  | 0 | <b>O</b> |  |
| One-color selection   |  | 0 | •        |  |

# Multiple Scan Patterns

The MC-500 Vixi has 22 preprogrammed scan patterns to allow treatment of varying retinal pathologies.



- \*1 For equal space patterns, No. v No. indicates the number of spots in horizontal and vertical directions.
- \*2 The arcade grid pattern is used for treatment of the periphery of macula in one-sixth units. The inner diameter is fixed and spot sizes range from 100 to 200 µm.

#### Typical Scan Patterns

Equal space (2v2, 3v3, 4v4, 5v5)

Square (2x2, 3x3, 4x4, 5x5)

The space between spots is equal in all directions.

The space between spots is equal in the horizontal and vertical directions.












The square pattern makes larger spaces in the diagonal direction than the horizontal and vertical directions.

The equal space pattern keeps the spaces between spots equal


allowing for denser photocoagulation than the square pattern.






# Auto Forward\*

Once photocoagulation is completed in one region, the MC-500 Vixi allows automated positioning of the scan pattern to the next region to undergo photocoagulation. This feature allows the surgeon to concentrate on focus adjustment.







The repeat mode with the auto forward function enables consecutive regions to undergo photocoagulation on a selected path without repeatedly pressing the foot switch.

<sup>\*</sup>The auto forward function is available for the equal space (2v2, 3v3, 4v4) and the square (2x2, 3x3, 4x4) patterns. The number of times of the forwarding differs depending on the scan pattern, spot size, and spacing.

# **Delivery Unit Options**

#### **Scan Delivery Units**



Scan slit lamp delivery unit (NIDEK SL-2000)



Scan slit lamp delivery unit (NIDEK SL-1800)



Scan attachable delivery unit (NIDEK SL-2000/SL-1800)



Scan attachable delivery unit (ZEISS SL 130)



Scan attachable delivery unit (HAAG 900BQ / 900BQ LED slit lamp)

#### **Single Delivery Units**



delivery unit (NIDEK SL-1800)



Attachable delivery unit (NIDEK SL-1800)



Attachable delivery unit (ZEISS SL 130)



BIO delivery unit (HEINE OMEGA 500)




YAG laser combination delivery unit (NIDEK YC-1800)




Endophotocoagulation delivery unit (ZEISS, LEICA)

#### **SOLIC (Safety Optics with Low Impact on Cornea)**


All scan slit lamp and slit lamp delivery units including attachable models incorporate the SOLIC optical design that ensures low energy density on the cornea and lens even for large spot sizes.




Small spot



Large spot





Large spot


#### Continuously Variable Spot Size

The spot size is continuously variable. The continuous variability enables the surgeon to easily compensate for the spot size change due to the use of a laser contact lens.



#### Stable Laser Power Output

Momentary increase followed by a plateau and an immediate decrease enables rapid and high-power laser emission on the scan patterns.



Time

# Practical and User-friendly Features

The intuitive graphic user interface and the easy-to-read color LCD touch screen allow easy and quick setup and confirmation of the scan pattern and treatment parameters.



The pop-up window appears once the displayed value, such as POWER, TIME or ITVL, is touched. It enables the surgeon to make significant changes to the laser values quickly with two-touch operations.

- 1. Touch the value on the screen
- 2. Select the value on the pop-up window



#### Memory of Scan Pattern

Four frequently used scan patterns can be saved and recalled with one-touch selection.

#### **Spot Size**

The spot size of the scan slit lamp delivery unit and slit lamp delivery unit is displayed on the LCD, and can be read with other parameters even in a dark room

#### **Actual Spot Size**

The converted spot size is displayed once the laser spot magnification of laser contact lens is selected on the pop-up window.



This screen displays all buttons and items for <mark>convenience s</mark>ake, but the actual screen is not consistent with this sample screen.

#### ■ Memory of Photocoagulation Data

In accordance with various clinical cases, up to 10 sets of photocoagulation data (color, power output, emission time, interval time, scan pattern, and spacing) can be registered.

Each set is retrievable quickly with one-touch operation.

| Memory No | Name |   | Color | Power | Time  | Intvl | Ptn   | SP   |
|-----------|------|---|-------|-------|-------|-------|-------|------|
|           |      |   | COIOI | mW    | sec   | sec   | Fui   | 31   |
| 1 PRP     | Scan | G | G     | 300   | 0.020 | 0.00  | 3 × 3 | 0.75 |
| 2 PRP     | Scan | Υ | Υ     | 300   | 0.020 | 0.00  | 3 v 3 | 0.75 |
| 3 P R P   | G    |   | G     | 200   | 0.200 | 0.40  |       |      |
| 4 PRP     | Υ    |   | Υ     | 150   | 0.200 | 0.40  |       |      |
| 5 BRV     | 0    |   | G     | 200   | 0.200 | 0.00  |       |      |
| 6 CRV     | О    |   | R     | 200   | 0.200 | 0.00  |       |      |
| 7 LI-1    |      |   | G     | 200   | 0.200 | 0.00  |       |      |
| 8 L I – 2 |      |   | G     | 1000  | 0.020 | 0.00  |       |      |
| 9         |      |   |       | 000   | 0.000 | 0.00  |       |      |
| 10        |      |   |       | 000   | 0.000 | 0.00  |       |      |
|           |      |   |       |       |       | ť     |       | Exit |


# LPM (Low Power Mode)

#### Minimally Invasive Photocoagulation

LPM (Low Power Mode) is a form of laser treatment that delivers reduced power to the retina. In LPM, the standard laser treatment power is reduced by a specified ratio.

To use the optional LPM, the MC-500 Vixi has to be equipped with the yellow laser. An additional software upgrade is required.

#### Laser Treatment for Macular Edema




Before treatment


After treatment with LPM Sample image – Outcomes may vary based on patient-specific response.

### Advantage of LPM Function

#### **Arcade Grid Scan Pattern**


In addition to regular mode, LPM includes a scan pattern that prevents treatment in a central circular area within the grid. Selecting the "Pos" (Position) button to align the aiming beam to the foveal center activates treatment that follows the grid scan pattern.





#### **Auto Forward Function**

The auto forward function is available for LPM. The MC-500 Vixi allows automated positioning of the scan pattern for photocoagulation.






#### **Main Body Specifications**

| 6 522  |
|--|
| Green: 532 nm                                |
| Yellow: 577 nm                               |
| Red : 647 nm                                 |
| Green: 50 to 1,700 mW*1                      |
| Yellow: 50 to 1,500 mW                       |
| Red: 50 to 800 mW*2                          |
| Continuous wave                              |
| 0.01 to 1.00 s, 2.00 s, 3.00 s*3             |
| 0.05 to 1.00 s*4                             |
| Laser diode, 670 nm, max. 0.6 mW             |
| Air cooling                                  |
| 115/230 V AC, 50/60 Hz                       |
| 400 VA                                       |
| 300 (W) x 480 (D) x 670 (H) mm / 35 kg       |
| 11.8 (W) x 18.9 (D) x 26.4 (H) " / 77.1 lbs. |
| 3D mouse, Control box, Remote control        |
|  |

- \*1 50 to 1,500 mW with scan delivery unit
- \*2 With the slit lamp delivery unit, scan slit lamp delivery unit, attachable delivery unit, and scan attachable delivery unit, the maximum power output is limited according to the spot size. Spot size 50  $\mu m$  - 500 mW, Spot size 60  $\mu m$  - 600 mW, Spot size 70  $\mu m$  - 700 mW
- \*3 0.01 to 0.05 second in scan mode
- \*4 0.10 to 1.00 second in auto manipulation mode and auto forward function



Foot print

#### Scan / Single Delivery Unit Specifications

| Model            | Scan delivery unit (MC-500 Vixi)  | Single delivery unit (MC-500)                          |
|------------------|---|--|
| Spot size        | 100 to 500 μm (scan mode, auto manipulation mode)                       | 50 to 1,000 μm (slit lamp, attachable deliveries)      |
|                  | 50 to 500 μm (single mode)  |  |
| Emission pattern | Single  | Single   |
|                  | Square (2x2, 3x3, 4x4, 5x5), Equal space (2v2, 3v3, 4v4, 5v5)*5,        |  |
|                  | Rectangle, Triple arc, Triple curve, Arcade grid, Triangle, Circle,     |  |
|                  | Arc (3/4 circle, 1/2 circle, 1/4 circle), Curve, Line                   |  |
| Туре             | Scan slit lamp delivery unit (NIDEK SL-2000/SL-1800)                    | Slit lamp delivery unit (NIDEK SL-1800)                |
|                  | Scan attachable delivery unit   | Attachable delivery unit (NIDEK SL-1800, ZEISS SL 130) |
|                  | (NIDEK SL-2000/SL-1800, ZEISS SL 130, HAAG 900BQ / 900BQ LED slit lamp) | BIO delivery unit (HEINE OMEGA 500)                    |
|                  |   | YAG laser combination delivery unit (NIDEK YC-1800)    |
|                  |   | Endophotocoagulation delivery unit (ZEISS, LEICA)      |
| Dimensions/mass  | 600 (W) x 450 (D) x 1,300 to 1,500 (H) mm / approximately 45 kg*6       | <b>←</b>   |
|                  | 23.6 (W) x 17.7 (D) x 51.2 to 59.1 (H)" / approximately 99.2 lbs.*6     | ←  |
|                  | (NIDEK SL-1800 scan slit lamp delivery with table)                      | (NIDEK SL-1800 slit lamp delivery with table)          |

- \*5 For equal space patterns, No. v No. indicates the number of spots in horizontal and vertical directions.
- \*6 The dimensions and mass differ depending on delivery types.





The laser wavelength of the purchased product is listed.



#### More clinical information available online at the NIDEK Education page

For more clinical information, please visit the Education page on the NIDEK website. This site allows access to case reports, journal articles, and video presentations.



Product/model name: Multicolor Laser Photocoagulator MC-500

Brochure and listed features of the device are intended for non-US practitioners.

Specifications may vary depending on circumstances in each country.

Specifications and design are subject to change without notice.

All brand and product names are trademarks or registered trademarks of their respective companies.



#### **HEAD OFFICE** (International Div.) 34-14 Maehama, Hiroishi-cho, Gamagori, Aichi 443-0038. JAPAN TEL: +81-533-67-8895 URL: www.nidek.com [Manufacturer]

#### TOKYO OFFICE (International Div.) 3F Sumitomo Fudosan Hongo

Bldg., 3-22-5 Hongo, Bunkyo-ku, Tokyo 113-0033, JAPAN TEL: +81-3-5844-2641 URL: www.nidek.com

#### NIDEK INC.

2040 Corporate Court, San Jose, CA 95131, U.S.A. TEL: +1-408-468-6400 +1-800-223-9044 (US Only)

URL: usa.nidek.com

#### NIDEK S.A.

Europarc, 13 rue Auguste Perret, 94042 Créteil, FRANCE TEL: +33-1-49 80 97 97 URL: www.nidek.fr

#### NIDEK TECHNOLOGIES S.R.L.

Via dell'Artigianato, 6/A, 35020 Albignasego (Padova), ITALY

TEL: +39 049 8629200 / 8626399 URL: www.nidektechnologies.it

#### NIDEK (SHANGHAI) CO., LTD. Rm3205,Shanghai Multi

Media Park, No.1027 Chang Ning Rd, Chang Ning District, Shanghai, CHINA 200050 TEL: +86 021-5212-7942 URL: www.nidek-china.cn

#### NIDEK SINGAPORE PTE. LTD.

51 Changi Business Park Central 2, #06-14, The Signature 486066, SINGAPORE TEL: +65 6588 0389 URL: www.nidek.sa